



The Ideal Skin Marker for withstanding alcoholic preparations

SurgiPEN

features a broad or precise fine line for marking during... General, Ortho, Neuro, Plastic, Micro and Vascular surgery, as well as radiology and radiation therapy.

SurgiPEN

- Surgical Marker, Gentian violet and Methylene blue dye, regular and fine tip, sterile, disposable, packed with 6 inch ruler.
- ► The marker's gentian violet and methylene blue ink is non-toxic, non-tattooing, and surgical prep-resistant.
- Marker's non-abrasive tip minimizes tissue trauma.
- A separate flexible ruler in the Marker kit is designed to conform to the patient's skin for easier, more precise marking.
- ▶ Ruler is calibrated in both inches and centimeters.
- Single use is recommended for best results.

Allows the surgeon to demarcate selected site, prior to or during surgical procedures. The skin is usually marked and left to dry for at least 60 seconds prior to scrubbing. Markings can be removed from skin with Isopropyl alcohol.

REGULAR TIP

Varicose Veins Thyroidectory Orthopaedics General Surgery Radiology etc.

FINE TIP

Plastic Surgery Micro Surgery Intricate general surgery facial, hand, etc., Paediatric surgery

ORDERING INFORMATION

Cat. No. Tip Ink Colour 25-12069 Regular Voilet 25-12070 Fine Voilet 25-12071 Regular Blue 25-12072 Fine Blue

Comes in packs of 25 sterile marker kits

Distributed by



For more information contact your Endosys Sales Representative +91-44-2826 9998, +91-93817 12055, +91-93845 60061 www.endosys.org | mail@endosys.org